


LIFE INSTRUMENTS®

McCULLOCH STYLE RETRACTOR



McCULLOCH RETRACTOR SYSTEM

Matte Finish Stainless Steel


744-2000-0

Standard Set

Catalog #	Quantity	Description
744-0001-0	1	Hinged Retractor Frame
749-2030-0	2	Small Muscle Blades 20x30mm
749-2040-0	2	Small Muscle Blades 20x40mm
749-2050-0	2	Small Muscle Blades 20x50mm
749-2060-0	2	Small Muscle Blades 20x60mm
749-2070-0	2	Small Muscle Blades 20x70mm
749-2080-0	2	Small Muscle Blades 20x80mm
749-2730-0	2	Large Muscle Blades 27x30mm
749-2740-0	2	Large Muscle Blades 27x40mm
749-2750-0	2	Large Muscle Blades 27x50mm
749-2760-0	2	Large Muscle Blades 27x60mm
749-2770-0	2	Large Muscle Blades 27x70mm
749-2780-0	2	Large Muscle Blades 27x80mm
749-0030-0	1	Hook 30mm
749-0040-0	1	Hook 40mm
749-0050-0	1	Hook 50mm
749-0060-0	1	Hook 60mm
749-0070-0	1	Hook 70mm
749-0080-0	1	Hook 80mm
744-2000-C	1	Sterilization Case for Retractor System

See reverse side for optional instruments

Components




FRAME

744-0001-0




HOOK



SMALL AND LARGE MUSCLE BLADES

Additional Items



**LONG SMALL
MUSCLE BLADE**

**LONG LARGE
MUSCLE BLADE**



STEPPED HOOK



LONG HOOK

Catalog #	Quantity	Description
749-2090-0	1	Long Small Muscle Blade 20x90mm
749-2000-0	1	Long Small Muscle Blade 20x100mm
749-2011-0	1	Long Small Muscle Blade 20x110mm
749-2012-0	1	Long Small Muscle Blade 20x120mm
749-2490-0	1	Long Large Muscle Blade 24x90mm
749-2400-0	1	Long Large Muscle Blade 24x100mm
749-2411-0	1	Long Large Muscle Blade 24x110mm
749-2412-0	1	Long Large Muscle Blade 24x120mm

Catalog #	Quantity	Description
749-0030-1	1	Stepped Hook 30mm
749-0040-1	1	Stepped Hook 40mm
749-0050-1	1	Stepped Hook 50mm
749-0060-1	1	Stepped Hook 60mm
749-0070-1	1	Stepped Hook 70mm
749-0080-1	1	Stepped Hook 80mm
749-0090-1	1	Stepped Hook 90mm
749-0100-1	1	Stepped Hook 100mm
749-0110-1	1	Stepped Hook 110mm
749-0090-0	1	Extra Long Hook 90mm
749-0100-0	1	Extra Long Hook 100mm
749-0110-0	1	Extra Long Hook 90mm