

DISK MARKERS



CERVICAL DISK MARKER

The markers easily identify the disk level on the x-ray.
Markers are 10cm in length with a 7mm protective ballstop.

Single	790-1070-0
Set of Three	790-1070-1

Designed by Edward Dunn, MD.



LUMBAR DISK MARKER

The markers easily identify the disk level on the x-ray.
Markers are 17.5cm in length with a 12.5mm protective ballstop.

Single	790-1013-0
Set of Three	790-1013-1

Saves Time • Reusable • Autoclavable

MARKERS



MARKERS

Round	790-1050-1
Square	790-1050-2

"WHITE" RUSSIAN FORCEP



Tamp End

"WHITE" RUSSIAN FORCEP

New forcep designed to facilitate collection, delivery, and packing of bone graft. Very useful for bone grafting adjacent to instrumentation, under rods and crosslinks.

Designed by Andrew White, MD.

8"	780-8005-0
10"	780-1005-0

JEFFERSON FORCEP



JEFFERSON FORCEP 7" length 175mm


1 x 2 Teeth	780-7012-0
Serrated	780-7013-0

LIFE INSTRUMENT ARM



**LIFE INSTRUMENT ARM WITH
UNIVERSAL ATTACHMENT (A)**
Table Clamp Included



400mm column 703-1601-0



**LIFE INSTRUMENT ARM WITH
LOW PROFILE ATTACHMENT (B)**
Holds Transverse Cervical Retractor Frame
Eliminates Cervical Retraction Movement During Operation
Table Clamp Included


703-1603-0

Transverse Retractor Frame not included.



702-0006-D Sterilization Insert to hold Arm
(703-1603-0 shown in tray above)

LIFE INSTRUMENT L-SHAPE COLUMN WITH ARTICULATING ARM



LIFE INSTRUMENT L-SHAPE COLUMN WITH ARTICULATING ARM Table Clamp Included

L-Shape Column (450x100mm) with Articulating Arm (400mm extension)

703-1801-0 with Universal Attachment

703-1803-0 with Low Profile Attachment*

* See page 10.3, figure B.

TABLE CLAMP



TABLE CLAMP

Stainless Steel

703-1600-CSS

FUNNEL



FUNNEL Stainless Steel

Bone Funnel Cup	706-0302-0
8.5" Tube	706-8508-1
6.5" Tube	706-6508-1

Cup disassembles from tube for easy cleaning.

* Sold Separately *

PLUNGER



PLUNGER

706-8508-2	for 8.5" Tube
706-6508-2	for 6.5" Tube

DURA & VESSEL REPAIR INSTRUMENTS



TITANIUM MICRO NEEDLEHOLDER

12" length 300mm

Curved with Lock Round Handle for Better Control

782-1201-1

8" length 200mm

782-0801-1



BAYONET MICRO NEEDLEHOLDER

120 working length

Curved with Lock

782-1201-3



ATRAUMATIC MICRO FORCEPS

Debaquey Style 1.2mm Round Handle for Better Control

12" length 300mm

780-1201-0

7" length 175mm

780-0801-0

** See page 11.21 for kit information.**

TO ORDER, CALL 800-925-2995

DURA & VESSEL REPAIR INSTRUMENTS



MICRO SCISSORS
Stainless Steel
Curved Round Handle for Better Control
12" length 300mm

783-1200-1

8" length 200mm

783-0800-1



BAYONET MICRO SCISSORS
Stainless Steel
Curved Round Handle for Better Control
120mm working length

783-1200-3

KNOT PUSHER



KNOT PUSHER

4" working length	784-1000-3
6" working length	784-1100-3
8" working length	784-1300-3



KNOT PUSHER

784-1105-3

Designed by Henry Ty, MD.

BAYONET BIPOLAR FORCEP



BAYONET BIPOLAR FORCEP 6" working length 1.0mm tip

780-1010-3	Straight Tip
780-1011-3	Angled Tip


KNIFE BLADE HOLDER



KNIFE BLADE HOLDER 6" working length

705-1010-3

MICRO SCISSOR




MICRO SCISSOR 7" length

801-0700-0

PERSONALIZED COBB ELEVATOR



All instruments can be personalized on handle with surgeon's name.




Great for Gifts!



Gift box included.

PATTIE TRAY



PATTIE TRAY

701-0507-0

DOTY-BOYD PATTIE TRAY



DOTY-BOYD PATTIE TRAY

701-0507-0-1

PATTIE PLATE



PATTIE PLATE

701-0305-0

PEDICLE PROBES



CURVED TAPERED TIP PROBE
Silicone Handle 10" overall length

525-0903-0



STRAIGHT TAPERED TIP PROBE
Silicone Handle 10" overall length


525-0913-0



PENFIELD STYLE PROBE
Silicone Handle 10" overall length

4mm	525-0904-1	slight curve
4mm	525-0904-0	flat

WIRE TWISTER



WIRE TWISTER

781-1000-1