

LIFE INSTRUMENTS®


RADIOLUCENT ZELPI® RETRACTOR



RADIOLUCENT ZELPI RETRACTOR

DEPTH	CATALOG #
2"	641-1220-2
3"	641-1230-2
4"	641-1240-2
5"	641-1250-2

- Removable Radiolucent Arms
- Excellent for Intraoperative Imaging
- Trademark spindle locking



RADIOLUCENT ZELPI RETRACTOR SET

See back for set components.

641-0000-0

*** NEW PRODUCT ADDITIONS ***



ZELPI HANDLE


641-1230-2-H



ZELPI ADAPTER ATTACHMENT
(Need two per retractor)

649-0001-0

Enables modularity of adding radiolucent blades of various length and width.



RADIOLUCENT ZELPI BLADE TRAY

641-0001-T



CARBON FIBER PEEK BLADES

SIZE	20mm width	25mm width
50mm	649-2050-0	649-2550-0
70mm	649-2070-0	649-2570-0
90mm	649-2090-0	649-2590-0
110mm	649-2111-0	649-2511-0

Additional blade sizes available upon request.



RADIOLUCENT ZELPI RETRACTOR SET

641-0000-0

Catalog #	Description	Depth	Quantity
641-1220-2-LA	Left Arm	2"	2
641-1220-2-RA	Right Arm	2"	2
641-1230-2-LA	Left Arm	3"	2
641-1230-2-RA	Right Arm	3"	2
641-1240-2-LA	Left Arm	4"	2
641-1240-2-RA	Right Arm	4"	2
641-1250-2-LA	Left Arm	5"	2
641-1250-2-RA	Right Arm	5"	2
641-0000-T	Tray		1
641-1230-2-H	Handle		2

Life Instruments' Radiolucent Zelpis are not visible on Xray and provide excellent retraction for Lumbar Spine procedures. They have been particularly helpful during A/P intraoperative imaging, allowing surgeons to see critical anatomy without interference. Avoiding manipulation, removal, and replacement of our retractors has allowed surgeons to experience less radiation exposure and shorter operative times.

The Radiolucent Zelpis also provide excellent exposure and retraction in Anterior Hip Approach surgeries.